# Using R and Bioconductor for proteomics data analysis

Laurent Gatto Computational Proteomics Unit

Atelier PROSPECTOM 19 Nov 2014, Grenoble, France

Version of this document: e8c452c [2014-11-19 12:55:42 +0000]

# Setup

The follow packages will be used throughout this documents. R version 3.1.1 or higher is required to install all the packages using BiocInstaller::biocLite.

```
library("mzID")
library("MSnID")
library("MSGFplus")
library("MSnbase")
library("rpx")
library("MLInterfaces")
library("pRoloc")
library("pRolocdata")
library("rTANDEM")
library("MSGFplus")
library("MSGFgui")
library("rols")
library("hpar")
```

The most convenient way to install all the tutorials requirement (and more related content), is to install RforProteomics with all its dependencies.

```
library("BiocInstaller")
biocLite("RforProteomics", dependencies = TRUE)
```

## Introduction

This tutorial illustrates R / Bioconductor infrastructure for proteomics. Topics covered focus on support for open community-driven formats for raw data and identification results, packages for peptide-spectrum matching, data processing and analysis:

- Exploring available infrastructure
- Mass spectrometry data
- Getting data from proteomics repositories
- Handling raw MS data
- Handling identification data
- MS/MS database search
- Analysing search results
- High-level data interface
- Quantitative proteomics
- Importing third-party quantitative data

- Data processing and analysis
- Statistical analysis
- Machine learning
- Annotation
- Other relevant packages/pipelines

Links to other packages and references are also documented. In particular, the vignettes included in the RforProteomics package also contains relevant material.

## Exploring available infrastructure

In Bioconductor version 3.0, there are respectively 65 proteomics, 44 mass spectrometry software packages and 7 mass spectrometry experiment packages. These respective packages can be extracted with the proteomicsPackages(), massSpectrometryPackages() and massSpectrometryDataPackages() and explored interactively.

```
library("RforProteomics")
pp <- proteomicsPackages()
display(pp)</pre>
```

## Mass spectrometry data

Type	Format	Package
raw	mzML, mzXML, netCDF, mzData	mzR (read)
identification	mzIdentML	${\tt mzR} \ {\rm and} \ {\tt mzID} \ ({\rm read})$
quantitation	mzQuantML	
peak lists	mgf	MSnbase (read/write)
other	mzTab	${\tt MSnbase} \ ({\rm read/write})$

## Getting data from proteomics repositories

Contemporary MS-based proteomics data is disseminated through the ProteomeXchange infrastructure, which centrally coordinates submission, storage and dissemination through multiple data repositories, such as the PRIDE data base at the EBI for MS/MS experiments, PASSEL at the ISB for SRM data and the MassIVE resource. The rpx is an interface to ProteomeXchange and provides a basic and unified access to PX data.

```
library("rpx")
pxannounced()
```

## ## 15 new ProteomeXchange annoucements

```
## Data.Set Publication.Data Message
## 1 PXD001283 2014-11-19 08:44:40 New
## 2 PXD001301 2014-11-19 08:42:01 New
## 3 PXD000837 2014-11-19 08:30:45 Updated information
## 4 PXD000715 2014-11-18 16:34:39 Updated information
```

```
## 5 PXD000837 2014-11-18 16:30:02 Updated information
## 6 PXD001354 2014-11-18 16:29:15 Updated information
## 7 PXD000627 2014-11-18 16:28:07 Updated information
## 8 PXD001125 2014-11-18 16:27:02 Updated information
## 9 PXD001045 2014-11-18 16:26:04 Updated information
## 10 PXD001260 2014-11-18 16:24:55 Updated information
## 11 PXD001414 2014-11-18 16:22:44 Updated information
## 12 PXD000715 2014-11-18 09:35:10
## 13 PXD000837 2014-11-18 09:27:36
## 14 PXD001260 2014-11-18 09:13:08 Updated information
## 15 PXD001045 2014-11-18 09:12:15 Updated information
px <- PXDataset("PXD000001")</pre>
px
## Object of class "PXDataset"
## Id: PXD000001 with 8 files
## [1] 'F063721.dat' ... [8] 'erwinia_carotovora.fasta'
## Use 'pxfiles(.)' to see all files.
pxfiles(px)
## [1] "F063721.dat"
## [2] "F063721.dat-mztab.txt"
## [3] "PRIDE_Exp_Complete_Ac_22134.xml.gz"
## [4] "PRIDE_Exp_mzData_Ac_22134.xml.gz"
## [5] "PXD000001 mztab.txt"
## [6] "TMT_Erwinia_1uLSike_Top10HCD_iso12_45stepped_60min_01.mzXML"
## [7] "TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.raw"
## [8] "erwinia_carotovora.fasta"
```

Other metadata for the px dataset:

```
pxtax(px)
pxurl(px)
pxref(px)
```

Data files can then be downloaded with the pxget function as illustrated below.

```
## mzf <- pxget(px, pxfiles(px)[6])
mzf <- "TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML"
mzf

## [1] "TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML"</pre>
```

#### Exercise

Explore what data files have been deposited by Pandey's recent draft map of the human proteome.

#### Solution

# Handling raw MS data

The mzR package provides an interface to the proteowizard code base, the legacy RAMP is a non-sequential parser and other C/C++ code to access various raw data files, such as mzML, mzXML, netCDF, and mzData. The data is accessed on-disk, i.e it does not get loaded entirely in memory by default. The three main functions are openMSfile to create a file handle to a raw data file, header to extract metadata about the spectra contained in the file and peaks to extract one or multiple spectra of interest. Other functions such as instrumentInfo, or runInfo can be used to gather general information about a run.

```
library("mzR")
ms <- openMSfile(mzf)</pre>
## Mass Spectrometry file handle.
## Filename: TMT_Erwinia_1uLSike_Top10HCD_iso12_45stepped_60min_01.mzXML
## Number of scans: 7534
hd <- header(ms)
dim(hd)
## [1] 7534
              21
names (hd)
##
    [1] "seqNum"
                                    "acquisitionNum"
##
    [3] "msLevel"
                                    "polarity"
##
    [5] "peaksCount"
                                    "totIonCurrent"
##
   [7] "retentionTime"
                                    "basePeakMZ"
   [9] "basePeakIntensity"
                                    "collisionEnergy"
## [11] "ionisationEnergy"
                                    "lowMZ"
  [13] "highMZ"
                                    "precursorScanNum"
## [15] "precursorMZ"
                                    "precursorCharge"
## [17] "precursorIntensity"
                                    "mergedScan"
## [19] "mergedResultScanNum"
                                    "mergedResultStartScanNum"
## [21] "mergedResultEndScanNum"
```

#### Exercise

Extract the index of the MS2 spectrum with the highest base peak intensity and plot its spectrum. Is the data centroided or in profile mode?

```
hd2 <- hd[hd$msLevel == 2, ]
i <- which.max(hd2$basePeakIntensity)
hd2[i, ]</pre>
```

#### Solution

```
5404
                          5404
## 5404
                                     2
                                                        275
                                                                2283283712
                                               1
        retentionTime basePeakMZ basePeakIntensity collisionEnergy
##
                                           354288224
## 5404
              2751.31 859.5032
                             lowMZ highMZ precursorScanNum precursorMZ
##
        ionisationEnergy
## 5404
                        0 100.5031 1995.63
                                                        5403
##
        precursorCharge precursorIntensity mergedScan mergedResultScanNum
## 5404
                                  627820480
##
        {\tt mergedResultStartScanNum\ mergedResultEndScanNum\ }
## 5404
pi <- peaks(ms, hd2[i, 1])</pre>
mz <- hd2[i, "basePeakMZ"]</pre>
par(mfrow = c(2, 2))
plot(pi, type = "h", main = paste("Acquisition", i))
plot(pi, type = "h", xlim = c(mz-0.5, mz+0.5))
pj <- peaks(ms, 100)
plot(pj, type = "l", main = paste("Acquisition", 100))
plot(pj, type = "l", xlim = c(536,540))
```

seqNum acquisitionNum msLevel polarity peaksCount totIonCurrent

Read the MSnbase::MSmap manual and look at the example to learn how the mzR raw data support can be exploited to generate maps of slides of raw MS data. (Note that the hd variable containing the raw data header was missing in MSnbase version < 1.14.1.)

# Handling identification data

##


The RforProteomics package distributes a small identification result file (see ?TMT\_Erwinia\_1uLSike\_Top10HCD\_isol2\_45ste that we load and parse using infrastructure from the mzID package.

```
## ## Number of scans: 5287 ## Number of PSM's: 5563
```

Various data can be extracted from the mzID object, using one the accessor functions such as database, scans, peptides, ... The object can also be converted into a data.frame using the flatten function.

#### Exercise

Is there a relation between the length of a protein and the number of identified peptides, conditioned by the (average) e-value of the identifications?



#### Solution

The mzR package also support fast parsing of mzIdentML files with the openIDfile function. Compare it, it terms of output and speed with mzID.

## reading TMT\_Erwinia.mzid.gz... DONE!

```
## user system elapsed
## 27.540 0.100 27.679
```

#### head(fid0)

```
spectrumid scan number(s) acquisitionnum passthreshold rank
##
## 1
     scan=5782
                           5782
                                           5782
                                                          TRUE
## 2
      scan=6037
                           6037
                                           6037
                                                          TRUE
                                                                   1
## 3
      scan=5235
                           5235
                                           5235
                                                          TRUE
                                                                   1
## 4
      scan=5397
                           5397
                                           5397
                                                          TRUE
                                                                   1
## 5
      scan=6075
                           6075
                                           6075
                                                          TRUE
## 6
      scan=5761
                                                          TRUE
                           5761
                                           5761
                                                                   1
     calculatedmasstocharge experimentalmasstocharge chargestate
## 1
                   1080.2321
                                             1080.2325
                                                                   3
## 2
                   1002.2115
                                             1002.2089
                                                                   3
## 3
                                                                   3
                   1189.2800
                                             1189.2836
## 4
                    960.5365
                                              960.5365
                                                                   3
                                                                   3
## 5
                   1264.3419
                                             1264.3409
                   1268.6501
                                             1268.6429
     ms-gf:denovoscore ms-gf:evalue ms-gf:rawscore ms-gf:specevalue
##
## 1
                    174 5.430080e-21
                                                 147
                                                          3.764831e-27
## 2
                    245 9.943751e-20
                                                          6.902626e-26
                                                 214
## 3
                    264 2.564787e-19
                                                 211
                                                          1.778789e-25
## 4
                    178 2.581753e-18
                                                  154
                                                          1.792541e-24
## 5
                    252 2.178423e-17
                                                  188
                                                          1.510364e-23
## 6
                    138 2.329453e-17
                                                  123
                                                          1.618941e-23
     assumeddissociationmethod isotopeerror isdecoy post pre end start
## 1
                            HCD
                                            0
                                                FALSE
                                                          S
                                                              R 84
                                                                        50
## 2
                            HCD
                                            0
                                                FALSE
                                                          R
                                                              K 315
                                                                       288
## 3
                            HCD
                                            0
                                                FALSE
                                                              R 224
                                                                       192
## 4
                            HCD
                                            0
                                                FALSE
                                                              R 290
                                                                       264
## 5
                            HCD
                                            0
                                                FALSE
                                                          F
                                                              R 153
                                                                       119
## 6
                                                              K 286
                            HCD
                                            \cap
                                                FALSE
                                                                       264
     accession length
                                                               description
       ECA1932
                   155
## 1
                                               outer membrane lipoprotein
## 2
       ECA1147
                   434
                                                            trigger factor
       ECA0013
                   295
## 3
                                       ribose-binding periplasmic protein
                   290
       ECA1731
                                                                  flagellin
       ECA1443
                   298
## 5
                            UTP--glucose-1-phosphate uridylyltransferase
##
  6
       ECA1444
                   468 6-phosphogluconate dehydrogenase, decarboxylating
##
                                   pepseq modified modification
## 1 PVQIQAGEDSNVIGALGGAVLGGFLGNTIGGGSGR
                                              FALSE
                                                             <NA>
## 2
                                              FALSE
            TQVLDGLINANDIEVPVALIDGEIDVLR
                                                             <NA>
## 3
       TKGLNVMQNLLTAHPDVQAVFAQNDEMALGALR
                                              FALSE
                                                             <NA>
             SQILQQAGTSVLSQANQVPQTVLSLLR
                                              FALSE
                                                             <NA>
    PIIGDNPFVVVLPDVVLDESTADQTQENLALLISR
                                              FALSE
                                                             <NA>
## 5
##
  6
                 WTSQSSLDLGEPLSLITESVFAR
                                              FALSE
                                                             <NA>
##
                                                       spectrumFile
## 1 TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML
## 2 TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML
## 3 TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML
## 4 TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML
## 5 TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML
## 6 TMT_Erwinia_1uLSike_Top10HCD_isol2_45stepped_60min_01.mzXML
```

```
## databaseFile
## 1 erwinia_carotovora.fasta
## 2 erwinia_carotovora.fasta
## 3 erwinia_carotovora.fasta
## 4 erwinia_carotovora.fasta
## 5 erwinia_carotovora.fasta
## 6 erwinia_carotovora.fasta

## 6 erwinia_carotovora.fasta

system.time({
    id1 <- openIDfile(f)
    fid1 <- psms(id1)
})

## Error in (function (classes, fdef, mtable) : unable to find an inherited method for function 'psms' :
## Timing stopped at: 0.237 0.001 0.237</pre>
```

## Error in head(fid1): error in evaluating the argument 'x' in selecting a method for function 'head':

# MS/MS database search

head(fid1)

While searches are generally performed using third-party software independently of R or can be started from R using a system call, the rtandem package allows one to execute such searches using the X!Tandem engine. The shinyTANDEM provides a interactive interface to explore the search results.

```
library("rTANDEM")
?rtandem
library("shinyTANDEM")
?shinyTANDEM
```

Similarly, the MSGFplus package enables to perform a search using the MSGF+ engine, as illustrated below:

A graphical interface to perform the search the data and explore the results is also available:

```
library("MSGFgui")
MSGFgui()
```

#### Exercise

Search TMT\_Erwinia\_1uLSike\_Top10HCD\_isol2\_45stepped\_60min\_01.mzXML against the fasta file from PXD000001 using, for example, MSGFplus/MSGFgui.

#### Solution

1. Get the fasta database:

```
## fas <- pxget(px, pxfiles(px)[8])
fas <- "erwinia_carotovora.fasta"</pre>
```

2. One could run MSGF+ from the command-line directly from R:

```
msgf <- system.file(package = "MSGFplus", "MSGFPlus", "MSGFPlus.jar")
system(paste0("java -jar ", msgf))
cmd <- paste("java -jar", msgf, "-protocol 2 -inst 1 -s", mzf, "-d", fas)
cmd</pre>
```

## [1] "java -jar /home/lg390/R/x86\_64-unknown-linux-gnu-library/3.1/MSGFplus/MSGFPlus/MSGFPlus.jar -pr

```
system(cmd)
```

or, use MSGFplus:

or, through the graphical user interface:

```
library("MSGFgui")
MSGFgui()
```

## Analysing search results

The MSnID package can be used for post-search filtering of MS/MS identifications. One starts with the construction of an MSnID object that is populated with identification results that can be imported from a data.frame or from mzIdenML files.

```
library("MSnID")
msnid <- MSnID(".")</pre>
```

```
## Note, the anticipated/suggested columns in the
## peptide-to-spectrum matching results are:
## ------
## accession
## calculatedMassToCharge
## chargeState
## experimentalMassToCharge
## isDecoy
## peptide
## spectrumFile
## spectrumID
```

```
PSMresults <- read.delim(system.file("extdata", "human_brain.txt",
                                       package="MSnID"),
                          stringsAsFactors=FALSE)
psms(msnid) <- PSMresults</pre>
show(msnid)
## MSnID object
## Working directory: "."
## #Spectrum Files: 1
## #PSMs: 997 at 37 % FDR
## #peptides: 687 at 57 % FDR
## #accessions: 665 at 65 \% FDR
The package then enables to define, optimise and apply filtering based for example on missed cleavages,
identification scores, precursor mass errors, etc. and assess PSM, peptide and protein FDR levels.
msnid$msmsScore <- -log10(msnid$`MS.GF.SpecEValue`)</pre>
msnid$absParentMassErrorPPM <- abs(mass_measurement_error(msnid))</pre>
filt0bj <- MSnIDFilter(msnid)</pre>
filtObj$absParentMassErrorPPM <- list(comparison="<", threshold=5.0)
filtObj$msmsScore <- list(comparison=">", threshold=8.0)
show(filt0bj)
## MSnIDFilter object
## (absParentMassErrorPPM < 5) & (msmsScore > 8)
filtObj.grid <- optimize_filter(filtObj, msnid, fdr.max=0.01,</pre>
                                  method="Grid", level="peptide",
                                  n.iter=500)
show(filtObj.grid)
## MSnIDFilter object
## (absParentMassErrorPPM < 2.3) & (msmsScore > 7.8)
msnid <- apply_filter(msnid, filt0bj.grid)</pre>
show(msnid)
## MSnID object
## Working directory: "."
## #Spectrum Files: 1
## #PSMs: 346 at 0 % FDR
## #peptides: 160 at 0 % FDR
## #accessions: 132 at 0 % FDR
```

The resulting data can be exported to a data.frame or to a dedicated MSnSet data structure for quantitative MS data, described below, and further processed and analyses using appropriate statistical tests.

# High-level data interface

The above sections introduced low-level interfaces to raw and identification results. The MSnbase package provides abstractions for raw data through the MSnExp class and containers for quantification data via the MSnSet class. Both store

- 1. the actual assay data (spectra or quantitation matrix), accessed with spectra (or the [, [[ operators) or exprs;
- 2. sample metadata, accessed as a data.frame with pData;
- 3. feature metadata, accessed as a data.frame with fData.

The figure below give a schematics of an MSnSet instance and the relation between the assay data and the respective feature and sample metadata.

Another useful slot is processingData, accessed with processingData(.), that records all the processing that objects have undergone since their creation (see examples below).

The readMSData will parse the raw data, extract the MS2 spectra (by default) and construct an MS experiment file

(Note that while readMSData supports MS1 data, this is currently not convenient as all the data is read into memory.)

```
## [1] "dummyiTRAQ.mzXML"
```

```
msexp <- readMSData(rawFile, verbose = FALSE)
msexp</pre>
```

```
## Object of class "MSnExp"
## Object size in memory: 0.2 Mb
## - - - Spectra data - - -
## MS level(s): 2
## Number of MS1 acquisitions: 1
## Number of MSn scans: 5
## Number of precursor ions: 5
## 4 unique MZs
## Precursor MZ's: 437.8 - 716.34
## MSn M/Z range: 100 2016.66
## MSn retention times: 25:1 - 25:2 minutes
## - - - Processing information - - -
## Data loaded: Wed Nov 19 14:44:21 2014
  MSnbase version: 1.14.0
## - - - Meta data - - -
## phenoData
##
    rowNames: 1
##
     varLabels: sampleNames
     varMetadata: labelDescription
## Loaded from:
```

```
dummyiTRAQ.mzXML
## protocolData: none
## featureData
     featureNames: X1.1 X2.1 ... X5.1 (5 total)
##
     fvarLabels: spectrum
##
     fvarMetadata: labelDescription
## experimentData: use 'experimentData(object)'
MS2 spectra can be extracted as a list of Spectrum2 objects with the spectra accessor or with the [ operator.
Individual can be accessed with [[.
length(msexp)
## [1] 5
msexp[[2]]
## Object of class "Spectrum2"
## Precursor: 546.9586
## Retention time: 25:2
## Charge: 3
## MSn level: 2
## Peaks count: 1012
## Total ion count: 56758067
The identification results stemming from the same raw data file can then be used to add PSM matches.
fData(msexp)
##
        spectrum
## X1.1
## X2.1
               2
## X3.1
               3
## X4.1
## X5.1
## find path to a mzIdentML file
identFile <- dir(system.file(package = "MSnbase", dir = "extdata"),</pre>
                  full.name = TRUE, pattern = "dummyiTRAQ.mzid")
basename(identFile)
## [1] "dummyiTRAQ.mzid"
msexp <- addIdentificationData(msexp, identFile)</pre>
## reading dummyiTRAQ.mzid... DONE!
```


fData(msexp)

```
spectrum scan number(s) passthreshold rank calculatedmasstocharge
## X1.1
                                           TRUE
                                                    1
                                                                     645.0375
               1
                                1
## X2.1
               2
                                2
                                           TRUE
                                                                     546.9633
                                                    1
## X3.1
               3
                                             NA
                              NA
                                                   NΔ
                                                                           NA
## X4.1
                4
                               NA
                                             NA
                                                   NA
                                                                           NA
                                           TRUE
## X5.1
               5
                                5
                                                                     437.2997
                                                    1
        experimentalmasstocharge chargestate ms-gf:denovoscore ms-gf:evalue
## X1.1
                         645.3741
                                             3
                                                                77
                                                                       79.36958
## X2.1
                         546.9586
                                             3
                                                                39
                                                                       13.46615
## X3.1
                                NA
                                            NA
                                                                NA
                                                                              NA
## X4.1
                                NA
                                            NA
                                                                NA
                                                                              NA
                         437.8040
                                             2
                                                                 5
                                                                      366.38422
## X5.1
        ms-gf:rawscore ms-gf:specevalue assumeddissociationmethod
## X1.1
                            5.527468e-05
                    -39
                                                                  CID
## X2.1
                    -30
                            9.399048e-06
                                                                  CID
## X3.1
                     NA
                                       NA
                                                                 <NA>
## X4.1
                                       NA
                                                                 <NA>
                     NA
## X5.1
                    -42
                            2.577830e-04
                                                                  CID
                                                          accession length
        isotopeerror isdecoy post pre end start
## X1.1
                    1
                        FALSE
                                  Α
                                       R 186
                                                170 ECA0984; ECA3829
## X2.1
                    0
                        FALSE
                                  Α
                                       K
                                          62
                                                 50
                                                            ECA1028
                                                                        275
## X3.1
                 <NA>
                           NA <NA> <NA>
                                          NA
                                                 NA
                                                                <NA>
                                                                         NA
## X4.1
                 <NA>
                           NA <NA> <NA>
                                                                <NA>
                                                                         NA
                                          NA
                                                 NA
## X5.1
                        FALSE
                                  L
                                       K
                                          28
                                                 22
                                                            ECA0510
                                                                        166
##
                                                                             description
## X1.1 DNA mismatch repair protein; acetolactate synthase isozyme III large subunit
                  2,3,4,5-tetrahydropyridine-2,6-dicarboxylate N-succinyltransferase
## X3.1
                                                                                    <NA>
## X4.1
                                                                                    <NA>
## X5.1
                            putative capsular polysacharide biosynthesis transferase
                    pepseq modified modification
                                                                databaseFile
## X1.1 VESITARHGEVLQLRPK
                              FALSE
                                                NA erwinia_carotovora.fasta
## X2.1
            IDGQWVTHQWLKK
                               FALSE
                                                NA erwinia_carotovora.fasta
## X3.1
                                  NA
                      <NA>
                                                NA
                                                                        <NA>
## X4.1
                      <NA>
                                  NA
                                                                        <NA>
## X5.1
                   LVILLFR
                              FALSE
                                                NA erwinia carotovora.fasta
##
        identFile nprot npep.prot npsm.prot npsm.pep
## X1.1
                 2
                       2
                                  1
                                            1
                                                      1
## X2.1
                 2
                       1
                                  1
                                            1
                                                      1
## X3.1
               NA
                      NA
                                 NA
                                           NA
                                                     NA
## X4.1
               NA
                      NA
                                 NA
                                           NA
                                                     NA
## X5.1
                 2
                       1
                                  1
                                            1
                                                      1
```

Spectra and (parts of) experiments can be extraced and plotted.

## msexp[[1]]


```
## Object of class "Spectrum2"
## Precursor: 645.3741
## Retention time: 25:1
## Charge: 3
## MSn level: 2
## Peaks count: 2921
## Total ion count: 668170086
```



# msexp[1:3]

```
## Object of class "MSnExp"
## Object size in memory: 0.13 Mb
## - - - Spectra data - - -
## MS level(s): 2
## Number of MS1 acquisitions: 1
## Number of MSn scans: 3
## Number of precursor ions: 3
## 2 unique MZs
## Precursor MZ's: 546.96 - 645.37
```

```
## MSn M/Z range: 100 2016.66
## MSn retention times: 25:1 - 25:2 minutes
## - - - Processing information - - -
## Data loaded: Wed Nov 19 14:44:21 2014
## Data [numerically] subsetted 3 spectra: Wed Nov 19 14:44:21 2014
## MSnbase version: 1.14.0
## - - - Meta data - - -
## phenoData
##
    rowNames: 1
##
    varLabels: sampleNames
    varMetadata: labelDescription
## Loaded from:
    dummyiTRAQ.mzXML,
                         dummyiTRAQ.mzid
## protocolData: none
## featureData
     featureNames: X1.1 X2.1 X3.1
##
##
     fvarLabels: spectrum scan number(s) ... npsm.pep (30 total)
     fvarMetadata: labelDescription
##
## experimentData: use 'experimentData(object)'
plot(msexp[1:3], full=TRUE)
```



Coercion to a data.frame is straightforward.

```
as(msexp[[1]], "data.frame")[100:105, ]
```

```
## mz i
## 100 141.0990 588594.812
## 101 141.1015 845401.250
## 102 141.1041 791352.125
## 103 141.1066 477623.000
## 104 141.1091 155376.312
## 105 141.1117 4752.541
```

## Quantitative proteomics


There are a wide range of proteomics quantitation techniques that can broadly be classified as labelled vs. label-free, depending whether the features are labelled prior the MS acquisition and the MS level at which quantitation is inferred, namely MS1 or MS2.

	Label-free	Labelled
MS1	XIC	SILAC, 15N
MS2	Counting	iTRAQ, TMT

In terms of raw data quantitation, most efforts have been devoted to MS2-level quantitation. Label-free XIC quantitation has however been addressed in the frame of metabolomics data processing by the xcms infrastructure.

An MSnExp is converted to an MSnSet by the quantitation method. Below, we use the iTRAQ 4-plex isobaric tagging strategy (defined by the iTRAQ4 parameter; other tags are available).

```
plot(msexp[[1]], full=TRUE, reporters = iTRAQ4)
```



```
msset <- quantify(msexp, method = "trap", reporters = iTRAQ4, verbose=FALSE)
exprs(msset)</pre>
```

```
## X1.1 4483.320 4873.996 6743.441 4601.378

## X2.1 1918.082 1418.040 1117.601 1581.954

## X3.1 15210.979 15296.256 15592.760 16550.502

## X4.1 4133.103 5069.983 4724.845 4694.801

## X5.1 11947.881 13061.875 12809.491 12911.479
```

# processingData(msset)

```
## - - - Processing information - - -
```

```
## Data loaded: Wed Nov 19 14:44:21 2014
## iTRAQ4 quantification by trapezoidation: Wed Nov 19 14:44:23 2014
## MSnbase version: 1.14.0
```

Other MS2 quantitation methods available in quantify include the (normalised) spectral index SI and (normalised) spectral abundance factor SAF or simply a simple count method.

Note that spectra that have not been assigned any peptide (NA) or that match non-unique peptides (npsm > 1) are discarded in the counting process.

See also The isobar package supports quantitation from centroided mgf peak lists or its own tab-separated files that can be generated from Mascot and Phenyx vendor files.

Have a look at the **?quantify** documentation file and review the above by walking through the example.

## Importing third-party quantitative data

The PSI mzTab file format is aimed at providing a simpler (than XML formats) and more accessible file format to the wider community. It is composed of a key-value metadata section and peptide/protein/small molecule tabular sections.

```
## mztf <- pxget(px, pxfiles(px)[2])
mztf <- "F063721.dat-mztab.txt"
(mzt <- readMzTabData(mztf, what = "PEP"))

## Warning in readMzTabData(mztf, what = "PEP"): Support for mzTab version
## 0.9 only. Support will be added soon.

## Detected a metadata section
## Detected a peptide section

## MSnSet (storageMode: lockedEnvironment)
## assayData: 1528 features, 6 samples
## element names: exprs
## protocolData: none
## phenoData
## rowNames: sub[1] sub[2] ... sub[6] (6 total)</pre>
```

```
##
     varLabels: abundance
##
     varMetadata: labelDescription
## featureData
     featureNames: 1 2 ... 1528 (1528 total)
##
##
     fvarLabels: sequence accession ... uri (14 total)
     fvarMetadata: labelDescription
##
## experimentData: use 'experimentData(object)'
## Annotation:
## - - - Processing information - - -
## mzTab read: Wed Nov 19 14:44:23 2014
## MSnbase version: 1.14.0
```

It is also possible to import arbitrary spreadsheets as MSnSet objects into R with the readMSnSet2 function. The main 2 arguments of the function are (1) a text-based spreadsheet and (2) column names of indices that identify the quantitation data.

```
csv <- dir(system.file ("extdata" , package = "pRolocdata"),</pre>
           full.names = TRUE, pattern = "pr800866n_si_004-rep1.csv")
getEcols(csv, split = ",")
    [1] "\"Protein ID\""
                                       "\"FBgn\""
##
##
   [3] "\"Flybase Symbol\""
                                       "\"No. peptide IDs\""
    [5] "\"Mascot score\""
                                       "\"No. peptides quantified\""
##
##
   [7] "\"area 114\""
                                       "\"area 115\""
                                       "\"area 117\""
  [9] "\"area 116\""
##
## [11] "\"PLS-DA classification\""
                                       "\"Peptide sequence\""
## [13] "\"Precursor ion mass\""
                                       "\"Precursor ion charge\""
## [15] "\"pd.2013\""
                                       "\"pd.markers\""
ecols <- 7:10
res <- readMSnSet2(csv, ecols)
head(exprs(res))
##
     area.114 area.115 area.116 area.117
## 1 0.379000 0.281000 0.225000 0.114000
## 2 0.420000 0.209667 0.206111 0.163889
## 3 0.187333 0.167333 0.169667 0.476000
## 4 0.247500 0.253000 0.320000 0.179000
## 5 0.216000 0.183000 0.342000 0.259000
```

#### head(fData(res))

## 6 0.072000 0.212333 0.573000 0.142667

```
##
     Protein.ID
                        FBgn Flybase.Symbol No..peptide.IDs Mascot.score
## 1
        CG10060 FBgn0001104
                                G-ialpha65A
                                                                    179.86
                                                           3
                                                           5
## 2
        CG10067 FBgn0000044
                                     Act57B
                                                                    222.40
## 3
        CG10077 FBgn0035720
                                    CG10077
                                                           5
                                                                    219.65
## 4
        CG10079 FBgn0003731
                                                           2
                                       Egfr
                                                                     86.39
## 5
        CG10106 FBgn0029506
                                    Tsp42Ee
                                                           1
                                                                     52.10
        CG10130 FBgn0010638
## 6
                                  Sec61beta
                                                           2
                                                                     79.90
     No..peptides.quantified PLS.DA.classification Peptide.sequence
## 1
                                                  PM
```

```
## 2
                             9
                                                   PM
## 3
                             3
## 4
                             2
                                                   PM
                                                              GGVFDTIQK
## 5
                             1
## 6
                             3
                                             ER/Golgi
                                                    pd.2013 pd.markers
##
     Precursor.ion.mass Precursor.ion.charge
                                                                unknown
## 1
                                                         PM
                                                         PM
                                                                unknown
## 2
## 3
                                                    unknown
                                                                unknown
## 4
                                                                unknown
                                                         PM
                                              2 Phenotype 1
## 5
                 626.887
                                                                unknown
## 6
                                                                      ER
                                                   ER/Golgi
```

## Data processing and analysis

## Raw data processing


For raw data processing look at MSnbases's clean, smooth, pickPeaks, removePeaks and trimMz for MSnExp and spectra processing methods.

The MALDIquant and xcms packages also features a wide range of raw data processing methods on their own ad hoc data instance types.

#### Processing and normalisation

Each different types of quantitative data will require their own pre-processing and normalisation steps. Both isobar and MSnbase allow to correct for isobaric tag impurities normalise the quantitative data.

```
data(itraqdata)
qnt <- quantify(itraqdata, method = "trap",</pre>
                 reporters = iTRAQ4, verbose = FALSE)
impurities \leftarrow matrix(c(0.929,0.059,0.002,0.000,
                        0.020,0.923,0.056,0.001,
                        0.000,0.030,0.924,0.045,
                        0.000, 0.001, 0.040, 0.923),
                      nrow=4, byrow = TRUE)
## or, using makeImpuritiesMatrix()
## impurities <- makeImpuritiesMatrix(4)</pre>
qnt.crct <- purityCorrect(qnt, impurities)</pre>
processingData(qnt.crct)
## - - - Processing information - - -
## Data loaded: Wed May 11 18:54:39 2011
## iTRAQ4 quantification by trapezoidation: Wed Nov 19 14:44:25 2014
## Purity corrected: Wed Nov 19 14:44:25 2014
## MSnbase version: 1.1.22
plot0 <- function(x, y, main = "") {</pre>
    old.par <- par(no.readonly = TRUE)</pre>
    on.exit(par(old.par))
    par(mar = c(4, 4, 1, 1))
    par(mfrow = c(2, 2))
```



Various normalisation methods can be applied the MSnSet instances using the normalise method: variance stabilisation (vsn), quantile (quantiles), median or mean centring (center.media or center.mean), ...

```
qnt.crct.nrm <- normalise(qnt.crct, "quantiles")
plot0(qnt, qnt.crct.nrm)</pre>
```

The combineFeatures method combines spectra/peptides quantitation values into protein data. The grouping is defined by the groupBy parameter, which is generally taken from the feature metadata (protein accessions, for example).

```
## arbitraty grouping
g <- factor(c(rep(1, 25), rep(2, 15), rep(3, 15)))
prt <- combineFeatures(qnt.crct.nrm, groupBy = g, fun = "sum")</pre>
```

## Combined 55 features into 3 using sum

```
processingData(prt)
```

```
## - - - Processing information - - -
## Data loaded: Wed May 11 18:54:39 2011
## iTRAQ4 quantification by trapezoidation: Wed Nov 19 14:44:25 2014
## Purity corrected: Wed Nov 19 14:44:25 2014
## Normalised (quantiles): Wed Nov 19 14:44:26 2014
## Combined 55 features into 3 using sum: Wed Nov 19 14:44:26 2014
## MSnbase version: 1.1.22
```

Finally, proteomics data analysis is generally hampered by missing values. Missing data imputation is a sensitive operation whose success will be guided by many factors, such as degree and (non-)random nature of the missingness. Missing value in MSnSet instances can be filtered out and imputed using the filterNA and impute functions.

```
set.seed(1)
qnt0 <- qnt
exprs(qnt0)[sample(prod(dim(qnt0)), 10)] <- NA
table(is.na(qnt0))

##
## FALSE TRUE</pre>
```

```
qnt00 <- filterNA(qnt0)
dim(qnt00)</pre>
```

```
## [1] 44 4
```

209

```
qnt.imp <- impute(qnt0)
plot0(qnt, qnt.imp)</pre>
```




#### Exercise

The mzt instance created from the mzTab file has the following is a TMT 6-plex with the following design:

In this TMT 6-plex experiment, four exogenous proteins were spiked into an equimolar *Erwinia* carotovora lysate with varying proportions in each channel of quantitation; yeast enolase (ENO) at 10:5:2.5:12.5:10, bovine serum albumin (BSA) at 1:2.5:5:10:5:1, rabbit glycogen phosphory-lase (PHO) at 2:2:2:2:1:1 and bovin cytochrome C (CYT) at 1:1:1:1:1:2. Proteins were then digested, differentially labelled with TMT reagents, fractionated by reverse phase nanoflow UPLC (nanoACQUITY, Waters), and analysed on an LTQ Orbitrap Velos mass spectrometer (Thermo Scientic).

Explore the mzt data using some of the illustrated functions. The heatmap and MAplot (see

Maplot function), taken from the RforProteomics vignette, have been produced using the same data.



## Statistical analysis

R in general and Bioconductor in particular are well suited for the statistical analysis of data. Several packages provide dedicated resources for proteomics data:

- MSstats: A set of tools for statistical relative protein significance analysis in DDA, SRM and DIA experiments.
- msmsTest: Statistical tests for label-free LC-MS/MS data by spectral counts, to discover differentially expressed proteins between two biological conditions. Three tests are available: Poisson GLM regression, quasi-likelihood GLM regression, and the negative binomial of the edgeR package.

## library(msmsTests)

## Error in library(msmsTests): there is no package called 'msmsTests'


```
data(msms.dataset)
## Warning in data(msms.dataset): data set 'msms.dataset' not found
msms.dataset
## Error in eval(expr, envir, enclos): object 'msms.dataset' not found
e <- pp.msms.data(msms.dataset)
## Error in eval(expr, envir, enclos): could not find function "pp.msms.data"
## Error in eval(expr, envir, enclos): object 'e' not found
null.f <- "y~batch"</pre>
alt.f <- "y~treat+batch"</pre>
div <- apply(exprs(e),2,sum)</pre>
## Error in exprs(e): error in evaluating the argument 'object' in selecting a method for function 'exp
res <- msms.edgeR(e,alt.f,null.f,div=div,fnm="treat")</pre>
## Error in eval(expr, envir, enclos): could not find function "msms.edgeR"
head(res)
## MSnSet (storageMode: lockedEnvironment)
## assayData: 1 features, 4 samples
    element names: exprs
## protocolData: none
## phenoData: none
## featureData
##
     featureNames: 1
     fvarLabels: Protein.ID FBgn ... pd.markers (12 total)
    fvarMetadata: labelDescription
##
## experimentData: use 'experimentData(object)'
## Annotation:
## - - - Processing information - - -
## Subset [888,4][1,4] Wed Nov 19 14:44:26 2014
## MSnbase version: 1.14.0
```


• isobar also provides dedicated infrastructure for the statistical analysis of isobaric data.

## Machine learning

The MLInterfaces package provides a unified interface to a wide range of machine learning algorithms. Initially developed for microarray and ExpressionSet instances, the pRoloc package enables application of these algorithms to MSnSet data.


```
library("MLInterfaces")
library("pRoloc")
library("pRolocdata")
data(dunkley2006)
traininds <- which(fData(dunkley2006)$markers != "unknown")</pre>
ans <- MLearn(markers ~ ., data = t(dunkley2006), knnI(k = 5), traininds)
ans
## MLInterfaces classification output container
## The call was:
## MLearn(formula = markers ~ ., data = t(dunkley2006), .method = knnI(k = 5),
      trainInd = traininds)
## Predicted outcome distribution for test set:
##
##
       ER lumen ER membrane
                                      Golgi Mitochondrion
                                                               Plastid
##
             5
                          140
                                        67
                                                       51
                                                                    29
##
             PM
                     Ribosome
                                        TGN
                                                  vacuole
             89
                           31
                                         6
## Summary of scores on test set (use testScores() method for details):
     Min. 1st Qu. Median Mean 3rd Qu.
## 0.4000 1.0000 1.0000 0.9332 1.0000 1.0000
kcl <- MLearn( ~ ., data = dunkley2006, kmeansI, centers = 12)</pre>
## clusteringOutput: partition table
##
       2 3 4 5
##
                        6
                            7
                                8 9 10 11 12
                                9 112 44 31 157
## 56 29 55 49 29 78 40
## The call that created this object was:
## MLearn(formula = ~., data = dunkley2006, .method = kmeansI, centers = 12)
plot(kcl, exprs(dunkley2006))
```






Average silhouette width: 0.33


# PCA screeplot



# PCA colored by partition



```
hcl <- MLearn( ~ ., data = t(dunkley2006), hclustI(distFun = dist, cutParm = list(k = 4)))
hcl</pre>
```



A wide range of classification algorithms are also available, as described in the ?MLearn documentation page. The pRoloc package also uses MSnSet instances as input and ,while being conceived with the analysis of spatial/organelle proteomics data in mind, is applicable many use cases.

## Annotation

All the Bioconductor annotation infrastructure, such as biomaRt, GO.db, organism specific annotations, ... are directly relevant to the analysis of proteomics data. A total of 93 ontologies, including some proteomics-centred annotations such as the PSI Mass Spectrometry Ontology, Molecular Interaction (PSI MI 2.5) or Protein Modifications are available through the rols.

```
library("rols")
olsQuery("ESI", "MS")
```

```
## MS:1000073 MS:1000162
## "ESI" "HiRes ESI"
```

Data from the Human Protein Atlas is available via the hpar package.

## Other relevant packages/pipelines

- Analysis of post translational modification with isobar.
- Analysis of label-free data from a Synapt G2 (including ion mobility) with synapter.
- Analysis of spatial proteomics data with pRoloc.
- Analysis of MALDI data with the MALDIquant package.
- Access to the Proteomics Standard Initiative Common QUery InterfaCe with the PSICQUIC package.

Additional relevant packages are described in the RforProteomics vignettes.

#### Session information

```
## R version 3.1.1 Patched (2014-09-02 r66514)
## Platform: x86_64-unknown-linux-gnu (64-bit)
##
## attached base packages:
## [1] stats4
                 parallel stats
                                     graphics grDevices utils
                                                                    datasets
## [8] methods
                 base
##
## other attached packages:
## [1] lattice_0.20-29
                             hpar_1.8.0
                                                  rols_1.8.0
## [4] MSGFgui_1.0.1
                             rTANDEM_1.6.0
                                                  data.table_1.9.4
## [7] pRolocdata_1.5.2
                             pRoloc_1.7.1
                                                  MLInterfaces_1.46.0
## [10] cluster_1.15.3
                             annotate_1.44.0
                                                  XML_3.98-1.1
## [13] AnnotationDbi_1.28.1 GenomeInfoDb_1.2.3
                                                  IRanges_2.0.0
## [16] S4Vectors_0.4.0
                             rpx_1.2.0
                                                  MSGFplus_1.0.3
## [19] MSnID_1.0.0
                             mzID 1.4.1
                                                  RforProteomics_1.5.2
## [22] MSnbase_1.14.0
                             BiocParallel_1.0.0
                                                  mzR_2.0.0
## [25] Rcpp 0.11.3
                             Biobase 2.26.0
                                                  BiocGenerics_0.12.1
## [28] BiocInstaller_1.16.1 knitr_1.8
## loaded via a namespace (and not attached):
## [1] affy_1.44.0
                                     affyio_1.34.0
## [3] base64enc_0.1-2
                                     BatchJobs_1.5
## [5] BBmisc_1.8
                                     biocViews_1.34.1
##
   [7] BradleyTerry2_1.0-5
                                     brew_1.0-6
## [9] brglm_0.5-9
                                     car_2.0-21
## [11] caret_6.0-37
                                     Category_2.32.0
## [13] checkmate_1.5.0
                                     chron_2.3-45
## [15] class_7.3-11
                                     codetools_0.2-9
                                     DBI_0.3.1
## [17] colorspace_1.2-4
## [19] digest_0.6.4
                                     doParallel_1.0.8
## [21] e1071_1.6-4
                                     evaluate_0.5.5
                                     FNN_1.1
## [23] fail_1.2
## [25] foreach_1.4.2
                                     formatR_1.0
                                     genefilter_1.48.1
## [27] gdata 2.13.3
## [29] ggplot2_1.0.0
                                     graph_1.44.0
```

```
## [31] grid_3.1.1
                                     gridSVG_1.4-0
## [33] GSEABase_1.28.0
                                     gtable_0.1.2
## [35] gtools_3.4.1
                                     htmltools_0.2.6
## [37] httpuv_1.3.2
                                     impute_1.40.0
## [39] interactiveDisplay_1.4.0
                                     interactiveDisplayBase_1.4.0
## [41] iterators 1.0.7
                                     kernlab 0.9-19
## [43] labeling 0.3
                                     limma 3.22.1
## [45] lme4_1.1-7
                                     lpSolve_5.6.10
## [47] MALDIquant_1.11
                                     MASS_7.3-35
## [49] Matrix_1.1-4
                                     mclust_4.4
## [51] mime_0.2
                                     minqa_1.2.4
## [53] munsell_0.4.2
                                     mvtnorm_1.0-0
## [55] nlme_3.1-118
                                     nloptr_1.0.4
## [57] nnet_7.3-8
                                     pcaMethods_1.56.0
## [59] pls_2.4-3
                                     plyr_1.8.1
## [61] preprocessCore_1.28.0
                                     proto_0.3-10
## [63] proxy_0.4-13
                                     R6_2.0.1
## [65] randomForest 4.6-10
                                     RBGL 1.42.0
## [67] R.cache_0.10.0
                                     RColorBrewer_1.0-5
## [69] RCurl_1.95-4.3
                                     rda 1.0.2-2
## [71] reshape2_1.4
                                     rJava_0.9-6
## [73] RJSONIO_1.3-0
                                     R.methodsS3_1.6.1
## [75] R.oo_1.18.0
                                     rpart_4.1-8
## [77] RSQLite_1.0.0
                                     RUnit 0.4.27
## [79] R.utils_1.34.0
                                     sampling_2.6
## [81] scales_0.2.4
                                     sendmailR 1.2-1
## [83] sfsmisc_1.0-26
                                     shiny_0.10.2.1
## [85] shinyFiles_0.4.0
                                     splines_3.1.1
## [87] SSOAP_0.8-0
                                     stringr_0.6.2
## [89] survival_2.37-7
                                     tools_3.1.1
## [91] vsn_3.34.0
                                     xlsx_0.5.7
## [93] xlsxjars_0.6.1
                                     XMLSchema_0.7-2
## [95] xtable_1.7-4
                                     zlibbioc_1.12.0
```